

# The world is relying on faster decisions

## SARS-CoV-2 Rapid Antigen Test



Testing ready to go when you need to know  
**Delivering reliable results at the point of care.**

SARS-CoV-2 infections continue to increase. With each day seeing spikes in communities around the world, it is becoming more vital than ever to relieve frontline staff - our heroes in hospitals and clinics. Especially in areas where lab testing is limited or not available, there is an urgent need for decentralized diagnostics that provide fast results ready to go when you need to know. To control the pandemic, healthcare professionals and patients need to know today whether an infection is present or not.


- High test quality /performance
- No analyzer needed
- Fast/direct results while patient is waiting
- Easy workflow



Reading and understanding results of the Roche SARS-CoV-2 Rapid Antigen Test is simple. The control line (C) determines whether the test is valid and a colored test line (T) indicates the detection of antigens specific to SARS-CoV-2. You can read the result within 15 minutes for accurate results.

### Together Improving Lives

#### Rapid Antigen Tests


Item No.	Description	UOM	Price
711-9327592190	SARS-CoV-2 Rapid Antigen Test	25/kit	\$289.00
711-9338322190	SARS-CoV-2 Control Test	10/box	\$129.99



Control, test is valid  C  
 Antigens specific to SARS-CoV-2 detected  T



#### Positive




In case of a positive\* result, a colored test line T appears in the lower section of the result window.

#### Negative



In case of a negative\*\* result, only the control line C appears.

#### Invalid



If no control line C is visible, the test is always invalid.

**Learn more.**  
 Contact your Medline Representative